

# CHROMA™

# 3D

## MASTER SHADES

FOR  ORIGIN® *Beyond™ Plus*

## "Bulletproof Approach to Staining Zirconia."

*Nadine Payne - Payne Oral Design*



3D-Master shade matching for full contour zirconia restorations has long been a challenge for dental laboratories. The new CHROMA 3D-Master Shade liquids precisely match the Origin *Beyond Plus*™ zirconia for a quick and consistent approach to staining. Simply dip the restoration into the color liquid for 1 minute after applying the Origin Incisal Enhancer on the incisal area.

- Eliminate the tedious & inconsistent conversion to 3D-Master from the Classic Vita.
- Optimized for *Beyond Plus*™ zirconia for maximum translucency effects
- Water based, Non-acidic color liquid
- 29 Shades

### Chroma 3D Master Shades

|                                 |              |                         |
|---------------------------------|--------------|-------------------------|
| Individual Master Shade Bottles | 50 mL Bottle | \$39 <sup>00</sup> Each |
|---------------------------------|--------------|-------------------------|

### Chroma 3D Master Set

|                                |               |                     |
|--------------------------------|---------------|---------------------|
| Complete Master Set: 29 Shades | 50 mL Bottles | \$752 <sup>00</sup> |
|--------------------------------|---------------|---------------------|

\*Not compatible with conventional 3YS tetragonal zirconia



# ORIGIN®

B&D DENTAL  
TECHNOLOGIES

800.255.2839 Toll Free  
OriginCADCAM.com