

# Artimax™

## Does your hinge give you all the motion you need?

Without causing too much stress on your hands?  
Or without having to lock or unlock the hinges?

**Allows you to keep rotating natural excursions without using any force while still keeping case rigid**

The built-in hinge allows stress of lateral & protrusive movements without any separation as adhesive brands do

Very stable built-in hinge eliminates glue and remains very solid throughout the whole process

Natural linear protrusive & retrusive movement like that of a human jaw allowing a maximum of 8-9 mm for anterior edge-to-edge movements

Plastic legs allow quick bite adjustment (Heat & reset within 1 min) with the Artimax heat gun, leaving no black marks or bubbles

**Gives you the accurate & stable centric occlusion like that of metal articulators even after numerous motions**

*Hinge never pops open while attempting to make wide lateral movements or trying to follow the wear marks during the excursion movements*

**Easily allows ceramists to follow the wear marks of the opposing teeth in the most natural way without even losing any stability**

Vertical stop rod offers maximum stability and occlusion accuracy for free-end cases

Artimax New Full arch system uses smaller diameter pins for better accuracy especially for small anterior teeth

Longer and more flexible neck than conventional hinges provides linear protrusive & retrusive movement as a normal function of the mouth and jaw

**Natural Motion! Rigidity! & Stability!**

Only with **Artimax™** can you have all three...

1-800-255-2839

 **B&D** Dental Corporation™  
WE EVOLVE WITH YOUR NEEDS